SANITISERS & DISINFECTANTS SUPER-Q (406)







Super-Q (406)

CONCENTRATED QAC DISINFECTANT



Super-Q (406) is a concentrated QAC disinfectant. Suitable for soak sanitising, foot baths, sanitising wash-downs, hard surface wiping and food factories. Effective dilution 1:100 parts of water.

- Alcohol-free sanitiser suitable for Halal food establishments
- Rinse-free formula, non-toxic, safe to use after dilution
- Concentrated for greater economy, dilute up to 100 times
- Suitable to be used with fogging machines.

BS EN 1040:2005 TEST RESULTS FOR SUPER-Q

TEST: STAPHYLOCOCCUS AUREUS (ATCC 6538) - GRAM POSITIVE

Contact Time (Dilution)	Initial Count	Post Application	Log Reduction	% Kill
5 Minutes (1:100)	120 million	< 10	> 7.08	> 99.99999%

TEST: PSEUDOMONAS AERUGINOSA (ATCC 15442) - GRAM NEGATIVE

Contact Time (Dilution)	Initial Count	Post Application	Log Reduction	% Kill
5 Minutes (1:100)	140 million	< 10	> 7.15	> 99.99999%



Suitable for Factory Foot Baths



Suitable for Fogging Machines

SUPER CONCENTRATED

Dilute up to 1:100

	CERTIFICAT	TEST C		TE	PSB Singar) COTE
	CERTIFICADO +	SuperSteam Asia Pac 3013 Bedok Industrial #04-2128 Singapore 489979 FOR				
	RTI	Product:	Concentrated QAC Dis	infectant		
	5	Brand/Model:	Super-Q			
	лкат •	Detail:	The test micro-organis Staphylococcus aureus Pseudomonas aerugin	(ATCC 6538)		
	СЕРТИФИКАТ	Specification:	the evaluation of basi	and antiseptics. Quantitive s bactericidal activity of che method and requirements (ph	mical disinfectants	
	•	Test Report:	7191106605-CHM15-0	8-LYP '		
	御川	Date of Test Report:	03 Feb 2015			
	認証証書			and found to comply with	the requirements	
	٠	of the above standar	d(s).			
	CERTIFICATE	DR LIN JIANHUA	~			
	CERTI	Vice President (Chemi TÜV SÜD PSB Pte Ltd				
	٠	Certificate No: D 03366	late of Original Issue: 02 Mar 2015	Date of Last Revision: -	Date of Expiry: 01 Mar 2017	
12/10/1	ERTIFIKAT	property of TÜV SÜD PSI to the terms and condition				
No. 000 Mai	ZER		sched purschant to ene verms an	N VAN VYDERDON		TUV®

SPECIFICATIONS

Size	Super-Q (406)
5 L	0880-Super-Q-5L